



Recombinant Human Carbonic Anhydrase III

Catalogue #: RP-1000

Source: E. coli

Purity: > 90% by scanning SDS-PAGE

Molecular Weight: ~ 35 kD by SDS-PAGE

Purification Method: Combination of ion exchange and affinity chromatography

Storage Condition: -20 °C

Advantages: Contains no human or animal products

Uses: Calibrators and controls for immunoassays, western blot standards

Sizes Available:
0.05 mg
0.1 mg
1.0 mg
5.0 mg
10.0 mg

Bulk sizes and pricing available

*For Research Use Only.
Not for Human In-Vivo Use.*

International Point of Care Inc.
135 The West Mall, Unit 9
Toronto, Ontario, M9C 1C2
Canada

P: 416-847-2777

F: 416-847-2744

E: customerservice@ipocdx.com